


# the *portable* travel mobility scooter











#### **FOLDABLE & PORTABLE**

No matter where you want to go, you can takeyour Luggie with you. When you get there just **unfold**, **hop on**, and **off you go**!

With **lightweight lithium batteries** you can travel **up to 20 kms** between charges. And the battery charger will work in any country.

The Luggie is ideal for those who want to **travel the** world or just down to the local shops.

And it's so **easy to fold up**: you can pack it into the boot and take it with you.



#### **EASY TO LOAD**

Unlike other foldups you can half fold the Luggie for ease of loading by one person, or fully fold when space is at a premium. At 23 kgs it's truly portable.



#### **EASY TO ADJUST**

The Luggie can be adjusted in many configurations to suit your comfort and ease of use.







Since I bought the Luggie, I have been travelling to so many different places, and more often, because the small light scooter has been so easy to use.

Maureen's first Luggie lasted 7 years and she's now on her second, upgraded Luggie



The Luggie is such a great scooter for getting around the mall. And it's the best for loading into the car boot.

Kath





#### **DRIVING**

Easily folds up to fit in the boot of the car to go shopping, on holidays, or off to a show!

#### **CRUISING**

Take the Luggie on a cruise and get around the decks with ease.

#### **FLYING**

Fly with the Luggie – in Australia or overseas. The lithium-ion batteries are IATA approved for taking on as hand luggage and your Luggie travels free!

#### **TREKKING**

With eight different models there's sure to be a Luggie to suit you, even when the road is rough.

### WHY LUGGIE?

- Easy to fold up
- Easy to operate
- Easy to transport
- Cool design
- Highly engineered
- Eight models







My wife and I very much enjoy travelling and exploring foreign parts. I purchased a Luggie Scooter three and a half years ago and have used it extensively both at home and overseas.

Scott Sanderson in the USA on his Luggie



John Thearle is a polio survivor from the 1940's. His wife Mary-Lou loads the Luggie into the Camry boot and they have 'walked' around Brisbane parks during the COVID-19 crisis.



#### **OPTIONS:**

- Underseat carry bag
- Larger 16.5 a/h Li-ion battery
- Armrests for Eco and Standard
- Sheepskin seat cover
- Anti-scratch loading mat
- Travel Case



#### **SPECIFICATIONS:**

	Eco	Standard	Elite	Elite Plus	Super	Super Plus	Elite Deluxe	Super Deluxe
Dimensions (L x W), cm	95 x 45.5	95 x 45.5	95 x 45.5	99 x 45.5	98.9 x 52.5	98.9 x 52.5	99 x 58	99 x 58
Length folded, cm	65 cm	65 cm	65 cm	65 cm	67 cm	67 cm	N/A	N/A
Seat (L x W), cm	39.8 x 37.9	39.8 x 37.9	39.8 x 37.9	39.8 x 37.9	39.2 x 43.2	39.2 x 43.2	39.8 x 37.9	39.2 x 43.2
Ground Clearance	2.5 cm	2.5 cm	2.5 cm	2.5 cm	2.7 cm	2.8 cm	2.5 cm	2.9 cm
Maximum Load	114 kg	114 kg	145 kg	145 kg	163 kg	163 kg	145 kg	163 kg
Range#	15-25 km	10-20 km	15-25 km	15-25 km	15-23 km	20-30 km	15-25 km	20-30 km
Speed, Max forward	6 kph	6 kph	6 kph	6 kph	6.4 kph	8 kph	6.4 kph	8 kph
Weight excluding battery	23.5 kg	23.5 kg	23.5 kg	27.3 kg	26 kg	31.7 kg	29.5 kg	32.2 kg
Battery weight	12 kg	2 kg	2 kg	2 kg	2 kg	2 kg	2 kg	2 kg
Battery (amp hours)	12 Ah SLA	8.5 Ah Li-ion	10.5 Ah Li-ion	10.5 Ah Li-ion	10.5 Ah Li-ion	10.5 Ah Li-ion*	10.5 Ah Li-ion	10.5 Ah Li-ion*
Battery Charger	3 amp	3 amp	5 amp	5 amp	5 amp	5 amp	5 amp	5 amp
Tyres, front	152mm, solid	152mm, solid	152mm, solid	152mm, solid	152mm, solid	152mm, solid	178mm, solid	178mm, solid
Tyres, rear	178mm, solid	178mm, solid	178mm, solid	216mm, pneumatic	178mm, solid	228mm, pneumatic	216mm, pneumatic	228mm, pneumatic
Suspension, Front	X	X	Χ	<b>√</b>	Χ	<b>√</b>	√	<b>√</b>

## www.luggieaustralia.com.au



**IMPORTED BY:** Scooters Australia 4/355 South Gippsland Hwy Dandenong Vic 3175



<sup>#</sup> Range is an estimate only and is subject to terrain, payload and battery condition
\* 16.5 A/h battery optional extra
Maximum Load: Load capacity is subject to slope gradient. See website for more information.

Not to be used as a seat when being transported Specifications may change at any time without prior notice and should be regarded as a guide only